

#### **Bovine Albumin Peroxidase**

Catalog # ASR1056

# **Specification**

# **Bovine Albumin Peroxidase - Product Information**

Description BOVINE ALBUMIN Peroxidase conjugated

(BSA)

Conjugate Peroxidase (Horseradish)

Physical State
Host Isotype

Lyophilized
Albumin

Buffer 0.01 M Sodium Phosphate, 0.15 M Sodium

Chloride, pH 7.2

Species of Origin
Reconstitution Volume
Bovine
1.0 mL

Reconstitution Buffer Restore with deionized water (or

equivalent)

Stabilizer 10 mg/ml Polyethylene Glycol (PEG-8000)
Preservative 0.01% (w/v) Gentamicin Sulfate. Do NOT

add Sodium Azide!

## **Bovine Albumin Peroxidase - Additional Information**

# **Shipping Condition**

**Ambient** 

#### Purity

This product was prepared from normal serum by a multi-step process which includes delipidation, salt fractionation and selective followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase anti-Bovine Albumin and anti-Bovine Serum.

# **Storage Condition**

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

## **Precautions Note**

This product is for research use only and is not intended for therapeutic or diagnostic applications.

## **Bovine Albumin Peroxidase - Protein Information**

## **Bovine Albumin Peroxidase - Protocols**

Provided below are standard protocols that you may find useful for product applications.





- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

**Bovine Albumin Peroxidase - Images**